



SPECIALTY MEDICAL PRODUCTS, INC.

Billing: 50 Pennsylvania Ave Malvern, PA 19382

Shipping: 187 Pennsylvania Ave Malvern, PA 19382


Phone: (888) 698-1522 Fax (610) 644-3992

Email: info@smpcares.com

Website: www.smpcares.com

2600KL

Wheelchair ramp scale



2610KL

Wheelchair scale with dual ramps



Health o meter[®]
P r o f e s s i o n a l

Advantages

- ✓ **EMR Connectivity**
Scale interfaces with a computer, EMR software or other electronic device via various Health o meter[®] Professional Connectivity Solutions
- ✓ **Calculates Body Mass Index (BMI)**
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- ✓ **Versatile Platform**
Wide and low profile platform accommodates wheelchairs and also provides a comfortable weighing experience for bariatric patients
- ✓ **Interactive Display**
High-contrast color screen provides a user-friendly interface
- ✓ **Casters and Folding Ramps for Easy Portability and Storage**
4 wheels including 2 casters afford ease of movement and the folding ramps make it simple to store
- ✓ **On-Screen Help Menu**
Guides users through operating functions and provides easy access to multiple settings
- ✓ **Value-Added Features**
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option) included at no additional cost



Specifications

- Capacity: 1,000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 36" (w) x 32 1/4"(d) x 2 1/2" (h) / 914 mm (w) x 819 mm (d) x 64 mm (h)
- Ramp Size: 32 1/4" (w) x 9 7/8" (d) / 819 mm (w) x 248 mm (d)
- Display: 1 1/2" / 38 mm High-Contrast Color TFT-LCD Screen
- Wheels: 4 (2 casters)
- Power Source: 100-240V adapter included, 6 D batteries (not included)
- Connectivity: USB, optional Bluetooth®
- Optional Secondary Ramp for 2600KL, order item # B2600RAMP
- Functions: LB/KG Conversion, LB/KG Lock Out, BMI, Zero, Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off, Time/Date, Variable Auto Off Time, Audible/Mute Sound Option
- 2 Year Limited Warranty - Eligible for ScaleSurance Extended Warranty

Options & Accessories for the 2600KL & 2610KL

Model	Description
2600KL	Digital wheelchair ramp scale
2610KL	Digital wheelchair with dual ramps
SS-2600KL	ScaleSurance 2 year extended warranty
SS-2610KL	ScaleSurance 2 year extended warranty
ADPT30	Adapter included, item # for replacement or spare
B2600RAMP	Optional secondary ramp for 2600KL
BLUETOOTHONE-P	Wireless thermal printer

Healthometer® Professional



Advantages

- EMR Connectivity**
 Scale has the ability to interface with a computer, EMR software or other electronic device via various Healthometer® Professional Connectivity Solutions
- Motion Detection**
 Scale has the ability to weigh accurately despite movement of the subject on the scale
- Moveable Arm & Foot Rest**
 The flip-arm and foot rest allows for easier patient access and patient comfort
- Value-Added Features**
 Operating functions (LB/KG Conversion, LB/KG Lock Out, Zero Out/Tare, Auto Off, Reweigh) included at no additional cost



SPECIALTY MEDICAL PRODUCTS, INC.

Billing: 50 Pennsylvania Ave Malvern, PA 19382
Shipping: 187 Pennsylvania Ave Malvern, PA 19382
Phone: (888) 698-1522 Fax (610) 644-3992
Email: info@smpcares.com
Website: www.smpcares.com

Healthometer®
 Professional
www.homscales.com

See reverse side for additional specifications



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2595KL	32.625" x 24.875" x 41.875" / 19.67	100.25 lb	892076002626



Specifications

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4"(w) x 15 3/4" (d) x / 464 mm (w) x 400 mm (d)
- Display: 1 1/8" LCD
- Wheels: 2
- Connectivity: USB
- Power Source: 110-240V Adapter (included), 6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out /Tare, Auto Off



2595KL Display